

# TendonGel



## SOOTHES SWELLINGS AND BRUISES

FOR THE REGENERATION OF LIGAMENTS, TENDONS AND MUSCLES

Quick help for acute injuries  
to tendon and ligaments

Supports early resolution  
of bruising and edematous  
swelling

Aids the recovery of  
strained muscles



More information at: [www.derbymed.eu](http://www.derbymed.eu)



Sport



Recovery

## Product Information

TendonGel supports the self-healing of ligaments, tendons and muscles. The mixture of Arnica, Camphor, Menthol and Methyl soothes swelling after acute trauma. A cooling effect stimulates blood circulation and the regeneration of ligaments, tendons and muscles. The alcohol contained causes evaporation which is ideal for a quick cooling effect. The gel is easy to apply and goes through the skin and hair so it does not stick together.

### When and where to use:

- Areas with acute inflammation of the musculoskeletal system
- In case of injuries in ligaments and tendons
- As an aid to reduce localized swellings and bruises
  - Bruising (Hematoma)
  - Tendinitis
  - Tendovaginitis
  - Bursitis
- Muscle regeneration after intense exercise

## Application Process

Apply gel to the affected area until the swelling has subsided.


Apply up to 3 times per day.

Wrapping the affected area after massaging with the gel can enhance the soothing effect.

**TIP!**

Using derbymed Flexotendon at the same time optimally supports the tendons and ligaments.

# TendonGel



### Ingredients:

Aqua, Alcohol Denat., Isopropyl Alcohol, Propylene Glycol, Glycerin, Camphor, Menthol Carbomer, Mentha Arvensis Leaf Oil, PEG-40 Hydrogenated Castor Oil, Trideceth-9, Triethanolamine, Methyl Salicylate, Disodium Edta, Arnica Montana Flower Extract



Presentation: 250 ml

derbymed<sup>®</sup>  
is an exclusive brand of



**ARK**  
equine

Distributed IE & UK by:  
ARKequine  
Newbridge Industrial Estate  
Newbridge Co. Kildare, W12 YF63, Ireland  
Ph: 045 433552(IE) or 01292 800022 (UK)  
Email: enquiries@ark-equine.com  
www.ark-equine.com



**derbymed**<sup>®</sup>